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
Reported at 133. Kongress DGCH

The report on BESTA was presented at the 133. Congress of German Surgical Society in April 2016, Berlin

Reported at SESAM-2016

The report on BESTA was presented at the SESAM Conference, Lissabon, June 2016

BESTA Course as a Driver's License: Basic Endosurgical Simulation Training and Assessment

Maxim Gorshkov (1), Sergey Soyuzov (2), Nikolay Matveev (3), Lyubov Gribkova (1), Denis Gribkov (1)

1) Sechenov First Moscow State Medical University, Moscow, Russia
2) South Ural State Medical University, Chelyabinsk, Russia
3) Moscow State University of Medicine and Dentistry, Moscow, Russia

Background
Simulation methods of practical training have become an important part of the educational process in laparoscopic surgery. Attention to OR is less effective, attention of novices is scattered between the tasks that overtook them. Neither Russian Society of Endosurgeons (ROSED) nor Russian surgical, gynecological or urological societies have mandatory preclinical training and certification courses for basic endosurgical skills. Existing courses such as FLS, E-BLUS, LASTT, etc. do not meet all requirements of ROSMED and RSS. Our aim was to develop and validate a course for basic endosurgical training and assessment; identify principles of assessment and define credential criteria to allow residents to start their next stage of training in clinical environment, in OR.

Material and Methods
Russian Society for Simulation Education in Medicine ROSMED [rossomed] has established a Working group in April 2015. Experts in endosurgical simulation training (Gorshkov MD, Soyuzov SA, Matveev NL, Tsarkov PV, Shubina LB, Gribkov DM, Leoniev AB) participated in this panel. We analysed international standardized courses, such as: TopGun, Yates, MISTELS, FLS, PLUS, E-BLUS, SUIT, LASTT, as well as numerous programs for endoscopic surgery training conducted in Russia.
Routines endosurgical interventions in abdominal, urological and gynecological laparoscopy were segmented into 35 skills and manipulation. 15 skills were selected from this list for practical training. 10 of the simulation tasks were proposed not only for the training but as an assessment tool as well. Proficiency-based criteria were described for each certification exercise.
For validation of the selected tasks ROSMED has organized training and assessment sessions for endosurgeons and surgical residents participated in the XIX Convention of Russian Society of Endosurgeons (February 2016, Moscow) - see the Table.

Persons applied	20
Persons participated	19
Male	18
Female	1
Senior year students	10
Residents	9
Surgeons	0
Experienced	0
Experts in endosurgery	2
Some experience	3
No experience	13

Results
The Working Group of ROSMED formulated aims, structure, main characteristics and the list of skills of basic course of the simulation endosurgical training. Features: endosurgical, sic, practical orientation, teacher-replacement techniques, simulation, universal surgeons, urologists and gynaecologists, proficiency-based, without a fixed equipment, reproducible and standardized tasks, validity of training and assessment, objective and discriminant assessment.
Ten tasks have been selected for training and assessment and validated (see Table of the tasks at the right part of the poster). Both experienced surgeons and experienced residents performed the tasks of BESTA course, their results were collected and analysed. Face validity, content validity, construct validity and the validity of the course have been shown (see graphic below).

Conclusion
The Working Group of ROSMED have formulated principles of the universal prepping of basic endosurgical skills, its main characteristics, requirements and theoretical part; the list of training and certification tasks. 15 practical tasks and 15 practical exercises are combined in the course - Endosurgical Simulation Assessment (BESTA). The course is presented at the Congress of RSS and ROSMED and a training and assessment sessions should be proceeded further training in the clinical environment in OR after achieving the score of proficiency level only.

More info at: www.BESTA.guru

Scope Navigation
The task is to identify angles of the FLS system (20° angle of the upper handle and 30° angle of the lower handle) and identify 21 targets (table from the scope observation).

Flag Transfer
The task is to transfer the flag of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with both hands simultaneously.

Scope Enhancement/Correction
The task is to increase the angle of the upper handle and correct the angle of the lower handle. The task is to identify 21 targets (table from the scope observation) with the help of the observer.

Precision Cutting
The task is cutting the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Clip and Cut
The task is to identify the target of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Needle Guidance
The task is to identify the target of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Extracorporeal Suture
The task is to identify the target of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Establishing
The task is to identify the target of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Intercorporeal Knot Suture
The task is to identify the target of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Continuous Suture
The task is to identify the target of the FLS system (Fundamentals of Laparoscopic Surgery) from the upper handle to the lower handle with the help of the observer.

Reported at National Congress of Russian Surgical Society

Offered as a part of
State Accreditation of
surgical residents

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Reported at 134. Kongress DGCH

BESTA was discussed in details at the 134. Congress of German Surgical Society in March 2017, Munich

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Computer video analysis of the tasks of BESTA simulation training for automated proficiency assessment

Reported at SESAM-2017

The report on
BESTA Computer
Analysis was presented
at SESAM Conference,
Paris University,
Paris, June 2017

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


© The "Peg Transfer" is based on the Task 1 of the FLS Program. Fundamentals of Laparoscopic Surgery™ (FLS) Program is owned by Society of American Gastrointestinal and Endoscopic Surgeons (SAGES) and American College of Surgeons (ACS)

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The task is kind of endosurgical scavenger hunt. It is similar to the Task 1 (Scope Navigation), but in this time Trainee needs to use dissector in non-dominant hand to pull up the covers of each cell. That enables to visualize and recognize hidden from direct view numbers in the cells using the angled scope.

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
Aim of the task is to cut off a circle on the Training Gauze between 2 round marking lines. Cut should not touch the marking lines. The task should be practiced until it can be performed in 98 seconds with all cuts made within the section between the two lines of the Training Gauze; this level of performance should be achieved on 2 consecutive attempts.



© The Circle Cutting Task is based on the Task 2 of the FLS Program. Fundamentals of Laparoscopic Surgery™ (FLS) Program is owned by Society of American Gastrointestinal and Endoscopic Surgeons (SAGES) and American College of Surgeons (ACS)

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With the clip-applier in the dominant (right) hand clip is placed on the duplicated rubber around the proximal right pin of the rectangle. Use dissector for assistance. The next clip is placed around the far right pin of the hexagon. Then the two clips are placed on the rubbers of rectangle and hexagon connecting them together. Between both clips should be 10 mm or more. When applying the clips both jaws have to be visualized. Divide both rubbers by using Metzenbaum in two steps.

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

The red rubber band is stretched on the right pins in form of hexagon. Insert the needle with thread and stitch the rubber strip in the gap between right pins from outside inwards. Pull the thread through using the second needle-holder but leave the end of the thread in the ribbon. Repeat stitches moving counter-clockwise. When all six spaces between the pins will be sutured remove the instruments.

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Grasp the atraumatic suture by the thread and introduce into the trainer, place a long suture through the two marks in the penrose drain. Three single throws of a knot are tied extracorporeally. Using a knot pusher secure each throw onto the drain. At the end of the task, once all three throws have been secured onto the drain, cut both ends of the suture inside the trainer.


© The Task Extracorporeal Suture is based on the Task 4 of the FLS Program. Fundamentals of Laparoscopic Surgery™ Program is owned by Society of American Gastrointestinal and Endoscopic Surgeons (SAGES) and American College of Surgeons (ACS)

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Bring the pre-tied loop by the pusher to the foam and place it around the appendage at the mark. Secure the knot on the mark by sliding the pusher rod down. Once the knot is tied complete the task by cutting the end of loop suture material inside the trainer by Metzenbaum scissors.


© The Task Endo Loop is based on the Task 3 of the FLS Program. Fundamentals of Laparoscopic Surgery™ Program is owned by Society of American Gastrointestinal and Endoscopic Surgeons (SAGES) and American College of Surgeons (ACS)

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Introduce the atraumatic needle with 12 cm thread to the trainer and stitch it through the marks on the penrose drain. The first throw must be a surgeon's knot or double throw, followed by two single throws. Once all three throws have been secured onto the penrose drain, cut both ends of the suture inside the trainer.


© The Task 9. Intracorporeal Knot is based on the Task 5 of the FLS Program. Fundamentals of Laparoscopic Surgery™ Program is owned by Society of American Gastrointestinal and Endoscopic Surgeons (SAGES) and American College of Surgeons (ACS)

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Penrose drain has 6 markings. Introduce the atraumatic needle to the trainer and stitch it through the two marks on the penrose drain. Tie the first knot (surgeon's double throw, followed by two single throws). Passed needle through the rest 4 markings on a penrose drain. Secure suture by three throws at the end (2+1+1). Once second knot is tied, cut both ends of the suture inside the trainer.

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